



Arrow

Acute Haemodialysis Catheters & Large-Bore Catheters

Quality, Protection, Performance

Large-Bore Catheters with ARROWg^{ard} Antimicrobial Technology

Arrow Haemodialysis Catheters from Teleflex are equipped with the Blue FlexTip Catheter Tip that balances insertion capability with indwelling performance.

Arrow Uncoated and Coated Acute Haemodialysis You-Bend Catheters

| REF. | CATHETER | | | | SWG | | |
|-----------------------------------|----------|--------|----------------|------------------|--------|--------|-----|
| | LUMEN | Ø | LENGTH | INTERNAL LUMEN G | LENGTH | Ø | QTY |
| with You-Bend Large-Bore Catheter | | | | | | | |
| CU-13122-F | 2 | 12 Fr. | 5 1/8" (13 cm) | 12 - 12 | 60 cm | 0.035" | 5 |
| CU-23122-F | 2 | 12 Fr. | 5 1/8" (13 cm) | 12 - 12 | 60 cm | 0.035" | 5 |
| CV-13122-UF | 2 | 12 Fr. | 5 1/8" (13 cm) | 12 - 12 | 60 cm | 0.035" | 5 |
| with straight large-bore catheter | | | | | | | |
| CV-12122-F | 2 | 12 Fr. | 6" (16 cm) | 12 - 12 | 60 cm | 0.035" | 5 |
| CS-12122-F | 2 | 12 Fr. | 6" (16 cm) | 12 - 12 | 60 cm | 0.035" | 5 |
| CS-12123-F | 3 | 12 Fr. | 6" (16 cm) | 16 - 12 - 12 | 60 cm | 0.035" | 5 |
| CS-22123-F | 3 | 12 Fr. | 6" (16 cm) | 16 - 12 - 12 | 60 cm | 0.035" | 5 |

Arrow Acute Haemodialysis Catheter Facts and Information

Optional You-Bend
Extension Lines

ARROWg^{ard}
Antimicrobial Technology

Blue FlexTip
Catheter Tip

Rotating suture hub

Teleflex, the Teleflex logo, Arrow, ARROWg^{ard} and Blue FlexTip are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.
© 2017 Teleflex Incorporated. All rights reserved. 94 30 61 - AP PK 01 • REV A • 05 17 PDF