




## ARROW CVC SYSTEMS

Product List




# ARROW CVC SYSTEMS



With their ARROW brand, Teleflex Medical covers the complete range of central venous catheters. Arrow CVCs are available in an unmatched variety of lumen configurations, lengths, French and Gauge sizes, with either standard catheters or catheters with ARROWg<sup>+</sup>ard antimicrobial technology. Continuous development and optimisation of their products make ARROW a partner in whom you and your patients can place their trust.

## CONTENTS

Introduction.....	02
Single Lumen CVC.....	04
Double Lumen CVC .....	07
Triple Lumen CVC .....	12
Quad Lumen CVC .....	17
<b>Five Lumen CVC .....</b>	<b>22</b>
<b>European Arrow Select Kits.....</b>	<b>23</b>

 *latex-free: all catheters in this catalogue are free from latex*

## LUMEN DESIGNS



## SIZES

A 14/16	E 10/10	F 18/14	H 16/18/18	I 16/14/16
B 14/18	12/12	G 12/16/12	20/22/22	J 14/18/18/16
C 22/22	18/20		20/23/23	K 16/14/18/18/18
D 20/22				

# TELEFLEX MEDICAL BRINGS INNOVATIVE SOLUTIONS TO SHARP PROTECTION

## SHARPS AWAY LOCKING DISPOSAL CUP

It offers an easy, safe way of moving used sharps from procedure to sharps container. The small red cups secures and protects the needles stuck into it. With needle bevels safely protected, the Sharps Away cup can be disposed into your hospital's larger sharps container for final disposal.

## ARROW MAXIMAL BARRIER KITS ALSO FEATURE ARROWG<sup>+</sup>ARD BLUE PLUS<sup>®</sup>, MAXIMAL BARRIER AND SHARPS SAFETY COMPONENTS

The Maximal Barrier Insertion trays are designed to provide everything you need to ensure minimal risk and comply with the best practice guidelines from OSHA, INS, IHI, and the CDC. For the prevention of catheter-related infection by also including ARROWg<sup>+</sup>ard Blue PLUS catheter technology with antimicrobial impregnation both inside and out.

## LESS INTIMAL DAMAGE, IRRITATION AND REDUCED PERFORATION RISK


Most ARROW catheters feature the ArrowBlue FlexTip, an integral catheter tip which is more pliant than the catheter body. It remains patent, yet deflects in case of inadvertent contact with vessel wall.

## ARROW RAULERSON<sup>®</sup> INTRODUCER SYRINGE

The syringe has a hollow plunger/barrel containing a patented valving system, which allows a springwire guide to be placed directly through the syringe into the vessel lumen. This results in less trauma, less blood exposure risk and virtually no chance for air embolism.

## ARROW ADVANCER<sup>®</sup>

The ARROW Advancer saves time and money by helping you to easily straighten the J-tip of the spring-wire guide and insert it with only one hand, advancing into its proper position with your thumb.



# ARROWg<sup>+</sup>ard Blue

Over 30 studies and reviews prove ARROWg<sup>+</sup>ard catheter technology's ability to decrease the incidence of catheter-related bloodstream infection (CRBSI), reduce costs and save lives.<sup>1</sup>

For example, ARROWg<sup>+</sup>ard reduces bacterial colonisation of the catheter by 60 per cent and catheter-related infection by an even greater 80 per cent.<sup>2</sup>

And it has been proven effective against a wide array of gram-positive bacteria, including *S. epidermidis*, MRSA and *S. aureus*; gram-negative bacteria, including Enterococci and *Pseudomonas* strains; and fungi, including *C. albicans*.<sup>3</sup>

Our latest breakthrough, ARROWg<sup>+</sup>ard Blue PLUS, gives you that proven infection-fighting ability both outside and in. Both chlorhexidine and silver sulfadiazine continuously prep the subcutaneous catheter tract and kill organisms in places that you simply cannot reach.



## DRAMATIC COST SAVINGS

ARROWg<sup>+</sup>ard's proven infection-preventing power\* can save your institution hundreds of thousands of dollars per year.<sup>4</sup>


CVCS/YEAR	INFECTION COST WITH COMMODITY CVC	INFECTION COST WITH ARROWG+ARD BLUE**	TOTAL SAVINGS
500	\$ 547.620	\$ 140.024	\$ 407.596
1.000	\$ 1.095.240	\$ 280.048	\$ 815.192
2.000	\$ 2.190.480	\$ 560.096	\$ 1.630.384
5.000	\$ 5.476.200	\$ 1.400.240	\$ 4.075.960



**ARROWg<sup>+</sup>ard Blue® (A+g)**  
CVC antimicrobial impregnation with chlorhexidine/silver sulfadiazine on the outside of the catheter



**ARROWg<sup>+</sup>ard Blue PLUS® (AGB Plus)**  
CVC antimicrobial impregnation with chlorhexidine/silver sulfadiazine on the inside and outside of catheter, extension lines and hubs



**MOLDED, FLEXIBLE  
LOW-PROFILE  
JUNCTION HUBS**  
Designed for patient comfort  
and ease of maintenance

**FLEXIBLE  
POLYURETHANE**  
Designed to soften in  
situ and enhance indwelling  
characteristics

**BLUE FLEXTIP**  
Soft catheter tip balances  
insertion capability with  
indwelling performance

**ANTIMICROBIAL  
TECHNOLOGY**  
Patented treatment bonds  
chlorhexidine and silver  
sulfadiazine to the entire  
surface length



**PROXIMAL PORT**

- blood sampling
- medications
- blood administration



**MEDIAL PORT**

- total parenteral nutrition (TPN)
- medications (if TPN use is anticipated, use distal port instead)



**DISTAL PORT**

- CVP monitoring
- blood administration
- high-volume or viscous fluids
- colloids
- medication



**ADDITIONAL LUMEN(S)**

- infusion
- medication

# SEPARATE SINGLE LUMEN CVC

ARROW				
REF.	CATHETER	LUMEN	CE	QTY
ID-04200	16 Ga x 8"	A	X	25
ID-04301	16 Ga x 8"	A	X	25
ID-04700	14 Ga x 8"	A	X	25



LUMEN DESIGN



# 1 LUMEN CVC SETS AND KITS

## "ES" CVC SETS

### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire
- insertion needle 18 Ga x 6.35 cm
- dilator
- indwelling 18 Ga Cath. over 20 Ga needle
- syringe or Raulerson syringe

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
ES-04150	20 Ga x 5"	A	.025" J/SST		X	X	10
ES-04218	18 Ga x 8"	A	.025" J/SST		X	X	10
ES-04300	16 Ga x 8"	A	.032" J/SST		X	X	10
ES-04301	16 Ga x 8"	A	.032" J/SST	X	X	X	10
ES-04301P	16 Ga x 8"	A	.032" J/SST	X	X	X	10
ES-04306	16 Ga x 6"	A	.032" J/SST	X	X	X	10
ES-04400	16 Ga x 12"	A	.032" J/SST		X	X	10
ES-04522	22 Ga x 4"	A	.021" SST		X	X	10
ES-04700	14 Ga x 8"	A	.032" J/SST		X	X	10
ES-04701	14 Ga x 8"	A	.032" J/SST	X	X	X	10
ES-04701P	14 Ga x 8"	A	.032" J/SST	X	X	X	10
ES-04706	14 Ga x 6"	A	.032" J/SST	X	X	X	10
ES-04706P	14 Ga x 6"	A	.032" J/SST	X	X	X	10
ES-04730	14 Ga x 12"	A	.032" J/SST		X	X	10

## "CV" CVC SETS

### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire equipped by Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- syringe 5 ml

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-04706	14 Ga x 6"	A	.032" J/SST		X	X	10
CV-04701	14 Ga x 8"	A	.032" J/SST		X	X	10
CV-50014*	14 Ga x 8"	A	.032" J/SST		X	X	20
CV-50014BF**	14 Ga x 8"	A	.032" J/SST		X	X	20
CV-04306	16 Ga x 6"	A	.032" J/SST		X	X	10
CV-04301	16 Ga x 8"	A	.032" J/SST		X	X	10

\* without Blue FlexTip \*\* without extension line

**"CS" CVC SETS****SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire equipped by Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator, injection site cap
- indwelling 18 Ga Cath. over 20 Ga needle
- syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-.04225	16 Ga x 8"	A	.025" J/SST		X	X	10
CS-04300	16 Ga x 8"	A	.032" J/SST		X	X	10
CS-04400	16 Ga x 12"	A	.032" J/SST		X	X	10
CS-04700	14 Ga x 8"	A	.032" J/SST		X	X	10
CS-04701	14 Ga x 8"	A	.032" J/SST		X	X	10

**"AK" CVC SETS****SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle  
18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped with syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-04050	16 Ga x 16"	A	.032" SST			X	X	20
AK-04150-E	20 Ga x 5"	A	.025" J/SST			X	X	10
AK-04200	16 Ga x 8"	A	.032J/SST	St. Silk		X		10
AK-04210	16 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04211	16 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04212	16 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04225	16 Ga x 8"	A	.025" J/SST	St. Silk		X		10
AK-04250	16 Ga x 35"	A	.032" J/SST			X		10
AK-04300	16 Ga x 8"	A	.032J/SST	St. Silk		X		10
AK-04301	16 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04301-SP*	16 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04306	16 Ga x 6"	A	.032" J/SST	St.Nylon		X		10
AK-04306-SP*	16 Ga x 6"	A	.032" J/SST	St.Nylon		X		10
AK-04400	16 Ga x 12"	A	.032" J/SST	St. Silk		X		10
AK-04522	22 Ga x 4"	A	.021" J/SST	St. Silk		X		10
AK-04550	20 Ga x 5"	A	.025" J/SST	St. Silk		X		10
AK-04650-E	24Gax3.5"	A	.018" J/SST				X	10
AK-04660	24 Ga x 5"	A	.018" J/SST				X	10
AK-04700	14 Ga x 8"	A	.032" J/SST	St. Silk		X		10
AK-04800	16 Ga x 6"	A	.032" J/SST	St. Silk		X		10
AK-04825	16 Ga x 6"	A	.025" J/SST	St. Silk	X	X	X	10
CK-04711	14 Ga x 6"	A	.032" J/SST	St. Silk		X		10
SL-04710	14 Ga x 8"	A	.032" J/SST	St. Silk		X		10

\* Sharps safety components

**"CV" SINGLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SETS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- Raulerson Syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-24701-E	14 Ga x 8"	A	.032" J/SST	X	X	X	10
CV-24301-E	16 Ga x 8"	A	.032" J/SST	X	X	X	10

**"CS" SINGLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SETS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- indwelling 18 Ga Cath.  
over 20 Ga needle
- dilator
- Raulerson Syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-24301-E	16 Ga x 8"	A	.032" J/SST	X	X	X	10
CS-24306-E	16 Ga x 6"	A	.032" J/SST	X	X	X	10
CS-24701-E	14 Ga x 8"	A	.032" J/SST	X	X	X	10
CS-24706-E	14 Ga x 6"	A	.032" J/SST	X	X	X	10

**"AK" SINGLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SETS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- Raulerson Syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-24301	16 Ga x 8"	A	.032" J/SST	St. Silk	X	X		10
AK-24301-SK*	16 Ga x 8"	A	.032" J/SST	St. Silk	X	X		10
AK-24306	16 Ga x 6"	A	.032" J/SST	St. Silk		X		10


\* Sharps safety components



## SEPARATE DOUBLE LUMEN CVC

ARROW

REF.	CATHETER	LUMEN	CE	QTY
HO-14702	7 Fr x 16"	C	X	5
MC-12702	7 Fr x 6"	B	X	5
MC-16702	7 Fr x 8"	C	X	5
MC-17702	7 Fr x 8"	B	X	5

LUMEN  
DESIGN

## 2 LUMEN CVC SETS AND KITS

### "CV" DOUBLE LUMEN CVC SETS

#### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-12702	7 Fr x 6"	B	.025" J/SST		X	X	10
CV-16702	7 Fr x 8"	C	.032" J/SST		X	X	10
CV-17702	7 Fr x 8"	B	.032" J/SST		X	X	10
CV-16702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	10
CV-17702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	10
CV-14702	7 Fr x 12"	C	.032" J/SST		X	X	10
CV-15802	8 Fr x 8"	E	.032" J/SST		X	X	10

### "CS" DOUBLE LUMEN CVC SETS

#### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-12402	4 Fr x 2"	C	.018" J/SST		X	X	5
CS-12702-E	7 Fr x 6"	B	.025" J/SST	X	X	X	5
CS-12802	8 Fr x 6"	E	.032" J/SST		X	X	5
CS-14402	4 Fr x 5 1/8"	C	.018" J/SST		X	X	5
CS-14502	5 Fr x 5 1/8"	D	.021" J/SST		X	X	5
CS-15402-E	4 Fr x 3 1/8"	C	.018" J/SST		X	X	5
CS-15802-E	8 Fr x 8"	E	.032" J/SST	X	X	X	5
CS-16402	4 Fr x 12"	C	.018" J/SST		X	X	5
CS-16702	7 Fr x 8"	C	.032" J/SST		X	X	5
CS-16702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	5
CS-17702	7 Fr x 8"	B	.032" J/SST		X	X	5
CS-17702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	5

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-17752	7 Fr x 24"	B	.032" J/SST	X	X	X	5
CS-18702	7 Fr x 16"	B	.032" J/SST		X	X	5
ES-14402	4 Fr x 5 1/8"	C	.018" SST		X	X	5
ES-14702	7 Fr x 12"	C	.032" J/SST		X	X	5
HO-14702	7 Fr x 12"	C	.032" J/SST		X	X	5

## "AK" DOUBLE LUMEN CVC KITS

### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
AK-12402	4 Fr x 2"	C	.018" J/SST		X		5
AK-12702	7 Fr x 6"	B	.025" J/SST	X	X		5
AK-12702-SP*	7 Fr x 6"	B	.025" J/SST	X	X		5
AK-12702-CDC	7 Fr x 6"	B	.025" J/SST	X	X		5
AK-12802	8 Fr x 6"	E	.032" J/SST	X	X		5
AK-12802-SP*	8 Fr x 6"	E	.032" J/SST	X	X		5
AK-14412	4 Fr x 5 1/8"	C	.018" J/SST		X		5
AK-14502	5 Fr x 5 1/8"	D	.021" J/SST		X		5
AK-15402	4 Fr x 3 1/8"	C	.018" J/SST		X		5
AK-15802	8 Fr x 8"	E	.032" J/SST	X	X		5
AK-16402	4 Fr x 12"	C	.018" J/SST		X		5
AK-16702	7 Fr x 8"	C	.032" J/SST		X		5
AK-16702-CDC	7 Fr x 8"	C	.032" J/SST		X		5
AK-17702	7 Fr x 8"	B	.032" J/SST		X		5
AK-17702-CDC	7 Fr x 8"	B	.032" J/SST		X		5
AK-17702-SP*	7 Fr x 8"	B	.032" J/SST		X		5
BH-12702	7 Fr x 6"	B	.032" J/SST		X		10
CK-03000-B	7 Fr x 6"	C	.032" J/SST		X	X	5

\* Sharps safety components

## DOUBLE LUMEN CVC KITS: INSTRUMENT KITS

### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml
- drape, fenestrated drape
- suture, needle holder
- swabs
- gallipot

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
EU-12402-IMIN	4 Fr x 2"	C	.018" J/SST		X	X	5
EU-14502-IMIN	5 Fr x 5 1/8"	D	.021" J/SST	X	X	X	5
EU-15402-IMIN	4Fr x 3 1/8"	C	.018" J/SST	X	X	X	5
EU-12702-IMIN	7 Fr x 6"	B	.025" J/SST	X	X	X	5
EU-16702-IMIN	7 Fr x 8"	C	.032" J/SST	X	X	X	5
EU-12802-IMIN	8 Fr x 6"	E	.032" J/SST	X	X	X	5

\* Sharps safety components

**"CV" DOUBLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SETS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-22702-E	7 Fr x 5 1/8"	B	.025" J/SST	X	X	X	10
CV-27702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	10
CV-26702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	10
CV-25802-E	8 Fr x 8"	E	.032" J/SST	X	X	X	10

**"CS" DOUBLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SETS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-22502	5 Fr x 2"	D	.021" J/SST		X	X	5
CS-22702-E	7 Fr x 6"	B	.032" J/SST	X	X	X	5
CS-22802	8 Fr x 5 1/8"	E	.032" J/SST		X	X	5
CS-24402	4 Fr x 5 1/8"	C	.018" J/SST		X	X	5
CS-25502	5 Fr x 3 1/8"	D	.021" J/SST		X	X	5
CS-25802-E	8 Fr x 8"	E	.032" J/SST	X	X	X	5
CS-26702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	5
CS-27702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	5

**"CS" DOUBLE LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC SETS****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-46702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	5
CV-46702-E	7 Fr x 8"	C	.032" J/SST	X	X	X	5
CS-47702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	5
CV-47702-E	7 Fr x 8"	B	.032" J/SST	X	X	X	5
CS-45802-E	8 Fr x 8"	E	.032" J/SST	X	X	X	5
CS-42802	8 Fr x 6"	E	.032" J/SST		X	X	5
CV-45802-E	8 Fr x 8"	E	.032" J/SST		X	X	5

**"AK" DOUBLE LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC KITS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-22502	5 Fr x 2"	D	.021" J/SST	St. Silk		X		5
AK-22702	7 Fr x 6"	B	.025" J/SST	St. Silk	X	X		5
AK-22802	8 Fr x 6"	E	.032" J/SST	St. Silk	X	X		5
AK-22802-CDC*	8 Fr x 6"	E	.032" J/SST	St. Silk	X	X		5
AK-22802-SK*	8 Fr x 6"	E	.032" J/SST	St. Silk	X	X		5
AK-24402	4 Fr x 5 1/8"	C	.018" J/SST	St. Silk		X		5
AK-25502	5 Fr x 3 1/8"	D	.032" J/SST	St. Silk		X		5
AK-25802	8 Fr x 8"	E	.032" J/SST	St. Silk	X	X		5
AK-25802-SK*	8 Fr x 8"	E	.032" J/SST	St. Silk	X	X		5
AK-26702	7 Fr x 8"	C	.032" J/SST	St. Silk		X		5
AK-27702	7 Fr x 8"	B	.032" J/SST	St. Silk		X		5
NR-22802-B	8 Fr x 8"	E	.025" J/SST	St. Silk		X		5

\* Sharps safety components

**"AK" DOUBLE LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC KITS****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-42702	7 Fr x 6"	B	.025" J/SST	St. Silk	X	X		5
AK-42702-CDC (CHP)	7 Fr x 6"	B	.025" J/SST	St. Silk	X	X		5
AK-42702-SK*	7 Fr x 6"	B	.025" J/SST	St. Silk	X	X		5
AK-42802-SK	8 Fr x 5 1/8"	E	.025" J/SST	St. Silk	X	X		5
AK-45802-BSK (19 Ga Inw.)	8 Fr x 8"	E	.032" J/SST	St. Silk	X	X		5
AK-45802-SK (20 Ga Inw.)	8 Fr x 8"	E	.032" J/SST	St. Silk	X	X		5
AK-46702	7 Fr x 8"	C	.032" J/SST	St. Silk	X	X		5
AK-46702-SK*	7 Fr x 8"	C	.032" J/SST	St. Silk	X	X		5
AK-47702	7 Fr x 8"	B	.032" J/SST	St. Silk	X	X		5
AK-47702-CDC (CHP)	7 Fr x 8"	B	.032" J/SST	St. Silk	X	X		5
AK-47702-SK*	7 Fr x 8"	B	.032"	St. Silk	X	X		5

\* Sharps safety components

**MAXIMAL BARRIER DOUBLE LUMEN ARROW<sup>g</sup>ard BLUE PLUS (AGB+) CVC KITS****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- echogenic insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe and Raulerson 5 ml
- Maximal Sterile Barrier components, incl. full-body drape and Sharps Safety components

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-47702-MSB	7 Fr x 8"	B	0.032" J/SST	Curved Silk	X	X	X	5

**DOUBLE LUMEN ARROW<sup>g</sup>ard BLUE (A+g) CVC KITS: INSTRUMENT KITS****SET COMPONENTS**


- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-25502-IMIN	5 Fr x 3 1/8"	D	.032" J/SST	Curved Polyester		X	X	5

# SEPARATE TRIPLE LUMEN CVC

ARROW				
REF.	CATHETER	LUMEN	CE	QTY
ES-18703	7 Fr x 16"	H	X	5
HO-14703	7 Fr x 12"	H	X	5
MC-12123	12 Fr x 6"	G	X	5
MC-12703	7 Fr x 6"	H	X	5
MC-14703	7 Fr x 12"	H	X	5
MC-15123	12 Fr x 8"	G	X	5



LUMEN DESIGN



## 3 LUMEN CVC SETS AND KITS

### "CV" TRIPLE LUMEN CVC SETS

#### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-10753	7.5 Fr x 4"	I	0.032" J/SST		X	X	10
CV-12703	7 Fr x 6"	H	0.032" J/SST		X	X	10
CV-12703-E	7 Fr x 6"	H	0.032" J/SST	X	X	X	10
CV-12853	8.5 Fr x 6"	I	0.032" J/SST		X	X	10
CV-14703	7 Fr x 12"	H	0.032" J/SST		X	X	10
CV-14703-E	7 Fr x 12"	H	0.032" J/SST	X	X	X	10
CV-15703	7 Fr x 8"	H	0.032" J/SST		X	X	10
CV-15703-E	7 Fr x 8"	H	0.032" J/SST	X	X	X	10

### "CS" TRIPLE LUMEN CVC SETS

#### SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-12703	7 Fr x 6"	H	.032" J/SST		X	X	5
CS-12703-E	7 Fr x 6"	H	.032" J/SST	X	X	X	5
CS-12853	8.5 Fr x 6"	I	.032" J/SST		X	X	5
CS-14553	5.5 Fr x 12"	H	.018" J/SST		X	X	5
CS-14703	7 Fr x 12"	H	.032" J/SST		X	X	5
CS-14703-E	7 Fr x 12"	H	.032" J/SST	X	X	X	5
HO-14703	HandsOff	H	.032" J/SST		X	X	5
CS-14703-H	Heparin	H	.032" J/SST		X		5
CS-15553-E	5.5 Fr x 3 1/8"	H	.018" J/SST	X	X	X	5

**"BR" TRIPLE LUMEN STANDARD/HIGHFLOW CENTRAL VENOUS CATH. SETS W/ STOPCOCKS****SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
BR-10703	7 Fr x 4"	H	.032" J/SST		X	X	5
BR-10753	7.5 Fr x 4"	I	0.032" J/SST	X	X	X	5
BR-12703	7 Fr x 6"	H	.032" J/SST		X	X	5
BR-12853	8.5 Fr x 6"	I	.032" J/SST		X	X	5
BR-12853-E	8.5 Fr x 6"	I	.032" J/SST	X	X	X	5
BR-15853	8.5 Fr x 8"	I	.032" J/SST		X	X	5
BR-15853-E	8.5 Fr x 8"	I	.032" J/SST	X	X	X	5

**"BR" TRIPLE LUMEN STANDARD/HIGHFLOW ARROWg<sup>ard</sup> BLUE (A+g) CVC. SETS W/ STOPCOCKS****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml
- stopcocks

ARROW							
REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
BR-20753	7.5 Fr x 4"	I	0.032" J/SST		X	X	5
BR-24703	8.5 Fr x 12"	H	.032" J/SST		X	X	5

**"AK" TRIPLE LUMEN CVC KITS****SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-12703	7 Fr x 6"	H	.032" J/SST	St. Silk	X	X		5
AK-12703-B	7 Fr x 6"	H	.025" J/SST	St. Silk	X	X		5
AK-12703-CDC	7 Fr x 6"	H	.025" J/SST	St. Silk	X	X		5
AK-12703-SP*	7 Fr x 6"	H	.025" J/SST	St. Silk	X	X		5
AK-14553	5.5 Fr x 12"	H	.021" J/SST	St. Silk		X		5
AK-14703	7 Fr x 12"	H	.032" J/SST	St. Silk		X		5
AK-14703-A	7 Fr x 12"	H	.032" J/SST	St. Silk	X	X		5
AK-14703-H, heparin	7 Fr x 12"	H	.032" J/SST	St. Silk		X		5
AK-14703-SP*	7 Fr x 12"	H	.032" J/SST	St. Silk		X		5
AK-15553	5.5 Fr x 3 1/8"	H	.018" J/SST	St. Silk		X		5
AK-15703	7 Fr x 8"	H	.032" J/SST	St. Silk		X		5
AK-15703-A	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-CDC	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-ASP*	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-B	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-BCDC	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-BSP	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-CDC	7 Fr x 8"	H	.032" J/SST	St. Silk		X		5
AK-15703-H Heparin	7 Fr x 8"	H	.032" J/SST	St. Silk	X	X		5
AK-15703-SP*	7 Fr x 8"	H	.032" J/SST	St. Silk		X		5
AK-16553	5.5 Fr x 5 1/8"	H	.018" J/SST	St. Silk		X		5
NR-12703	7 Fr x 6"	H	.032" J/SST	St. Silk		X		5
YN-15703	7 Fr x 8"	H	.032" J/SST	St. Nylon		X		5

\* Sharps safety components

**"CV" TRIPLE LUMEN ARROWgard BLUE (A+g) CVC SET****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-20753	7.5 Fr x 4"	I	0.032" J/SST		X	X	10
CV-22703	7 Fr x 6"	H	0.032" J/SST		X	X	10
CV-22703-E	7 Fr x 6"	H	0.032" J/SST	X	X	X	10
CV-24703-E	7 Fr x 12"	H	0.032" J/SST	X	X	X	10
CV-25703	7 Fr x 8"	H	0.032" J/SST		X	X	10
CV-25703-E	7 Fr x 8"	H	0.032" J/SST	X	X	X	10



## TRIPLE LUMEN CVC: INSTRUMENT KITS

## SET COMPONENTS

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson® 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle
- drape, fenestrated drape
- suture, needle holder
- swabs
- gallipot

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-15553-IMIN	5,5 Fr x 3 1/8"	H	.018" J/SST	Curved Polyester		X	X	5
EU-12703-IMIN	7 FR x 6"	H	.032" J/SST	Curved Polyester	X	X	X	5
EU-15853-IMIN	7 FR x 8"	I	.032" J/SST	Curved Polyester	X	X	X	5
EU-12853-IMIN	8,5 Fr x 6"	I	.032" J/SST	Curved Polyester	X	X	X	5

"CS" TRIPLE LUMEN ARROWg<sup>ard</sup> BLUE (A+g) CVC SET

## SET COMPONENTS

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-22703-E	7 Fr x 6"	H	0.032" J/SST	X	X	X	5
CS-22853-E	8.5 Fr x 6"	I	0.032" J/SST	X	X	X	5
CS-24703-E	7 Fr x 12"	H	0.032" J/SST	X	X	X	5
CS-25703-E	7 Fr x 8"	H	0.032" J/SST	X	X	X	5
CS-25853	8.5 Fr x 8"	I	0.032" J/SST		X	X	5
CS-25853-E	8.5 Fr x 8"	I	0.032" J/SST	X	X	X	5

"AK" TRIPLE LUMEN ARROWg<sup>ard</sup> BLUE (A+g) CVC KITS

## SET COMPONENTS

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5ml

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-22703	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-24703	7 Fr x 12"	H	0.032" J/SST	St. Silk	X	X		5
AK-24703-A	7 Fr x 12"	H	0.032" J/SST	St. Silk	X	X		5
AK-24703-B	7 Fr x 12"	H	0.025" J/SST	St. Silk	X	X		5
AK-24703-ASK*	7 Fr x 12"	H	0.032" J/SST	St. Silk	X	X		5
AK-24703-SK*	7 Fr x 12"	H	0.032" J/SST	St. Silk	X	X		5
AK-25703	7 Fr x 8"	H	0.032" J/SST	St. Silk		X		5
AK-25703-A	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5
AK-25703-B	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5

\* Sharps safety components

**"AK" TRIPLE LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC KITS****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-42703	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-42703-C	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-42703-SK*	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-42703-AX*	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-42703-CDC	7 Fr x 6"	H	0.032" J/SST	St. Silk	X	X		5
AK-42703-BSK*	7 Fr x 6"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5
AK-45703-A	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5
AK-45703-ASK*	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5
AK-45703-AX*	7 Fr x 8"	H	0.032" J/SST	St. Silk	X	X		5
AK-45703-B	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-BSK*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-BX*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-C	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-AC	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-AVC	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-BC	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-CDC*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-ACDC*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-BCDC*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5
AK-45703-SK*	7 Fr x 8"	H	0.025" J/SST	St. Silk	X	X		5

\* Sharps safety components

**MAXIMAL BARRIER TRIPLE LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC KITS****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- echogenic insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ Raulerson syringe 5 ml
- Maximal Sterile Barrier components, incl. full-body drape and Sharps Safety components

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-45703-MSB	7 Fr x 8"	H	0.032" J/SST	Curved Silk	X	X	X	5

**“EU” TRIPLE LUMEN ARROWg<sup>ard</sup> BLUE (A+g) CVC KITS: INSTRUMENT KITS**

**SET COMPONENTS**


- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml
- drape, fenestrated drape
- suture, needle holder
- swabs
- gallipot

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-22703-IMIN	7Fr x 6"	H	.032" J/SST	Curved Polyester	X	X	X	5
EU-25703-IMIN	7 Fr x 6"	H	.032" J/SST	Curved Polyester	X	X	X	5

**SEPARATE QUAD LUMEN CVC**

LUMEN DESIGN



**4 LUMEN CVC SETS AND KITS**

**“CV” QUAD LUMEN CVC SET**

**SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-12854	8.5 Fr x 6"	J	.032" J/SST		X	X	10
CV-15854	8.5 Fr x 8"	J	.032" J/SST		X	X	10

**"CS" QUAD LUMEN CVC SET****SET COMPONENTS**

- polyurethane catheter (PUR) with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-12854-E	8.5 Fr x 6"	J	.032" J/SST	X	X	X	5
CS-14854-E	8.5 Fr x 12"	J	.032" J/SST	X	X	X	5
CS-15854-E	8.5 Fr x 8"	J	.032" J/SST	X	X	X	5

**"CV" QUAD LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SET****SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-22854	8.5 Fr x 6"	J	.032" J/SST		X	X	10
CV-25854	8.5 Fr x 8"	J	.032" J/SST		X	X	10

**"CV" QUAD LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC SET****SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle  
18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CV-42854	8.5 Fr x 6"	J	.032" J/SST		X	X	10
CV-45854	8.5 Fr x 8"	J	.032" J/SST		X	X	10

**“CS” QUAD LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC SET**

**SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-22854-E	8.5 Fr x 6"	J	.032" J/SST	X	X	X	5
CS-24854-E	8.5 Fr x 12"	J	.032" J/SST	X	X	X	5
CS-25854-E	8.5 Fr x 8"	J	.032" J/SST	X	X	X	5

**“CS” QUAD LUMEN ARROWg<sup>+</sup>ard BLUE PLUS (AGB+) CVC SET**

**SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-42854-E	8.5 Fr x 6"	J	.032" J/SST		X	X	5
CS-45854-E	8.5 Fr x 8"	J	.032" J/SST		X	X	5

**“AK” QUAD LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC KITS**

**SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-22854	8.5 Fr x 6"	J	.032" J/SST	St. Silk	X	X		5
AK-22854-SK*	8.5 Fr x 6"	J	.032" J/SST	St. Silk	X	X		5
AK-24854	8.5 Fr x 12"	J	.032" J/SST	St. Silk	X	X		5
AK-25854	8.5 Fr x 8"	J	.032" J/SST	St. Silk	X	X		5
AK-25854-SK*	8.5 Fr x 8"	J	.032" J/SST	St. Silk	X	X		5
NR-25854	8.5 Fr x 8"	J	.032" J/SST	St. Silk		X		5

\* Sharps safety components

**“AK” QUAD LUMEN ARROWg<sup>ard</sup> BLUE PLUS (AGB+) CVC KITS**

**SET COMPONENTS**

- polyurethane AGB + catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
AK-42854-CDC	8.5 Fr x 6"	J	.032" J/SST	St. Silk		X		5
AK-42854-SK*	8.5 Fr x 6"	J	.032" J/SST	St. Silk		X		5
AK-45854-CDC	8.5 Fr x 8"	J	.032" J/SST	St. Silk		X		5
AK-45854-SK*	8.5 Fr x 8"	J	.032" J/SST	St. Silk		X		5

\* Sharps safety components

**MAXIMAL BARRIER QUAD LUMEN ARROWg<sup>ard</sup> BLUE PLUS (AGB+) CVC KITS**

**SET COMPONENTS**

- polyurethane AGB+ catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- echogenic insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- Maximal Sterile Barrier components, incl. full-body drape and Sharps Safety components

ARROW

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-45854-MSB	8.5 Fr x 8"	J	.0032" J/SST	Curved Silk	X	X	X	5

**QUAD LUMEN ARROWg<sup>+</sup>ard BLUE (A+g) CVC KITS: INSTRUMENT KITS**

ARROW

**SET COMPONENTS**

- polyurethane A+g catheter with Blue FlexTip
- guide wire w/ or w/o Advancer
- insertion needle 18 Ga x 6.35 cm
- indwelling 18 Ga Cath. over 20 Ga needle
- Syringe or Raulerson 5 ml
- w/ or w/o local anesth., w/ equipped by syr + needle

REF.	CATHETER	LUMEN	SWG	SUTURE	RAULERSON	DILATOR	CE	QTY
EU-25854-IMIN	8.5 Fr x 6"	J	.032" J/SST	Curved Polyester	X	X	X	5

# SEPARATE FIVE LUMEN CVC

LUMEN DESIGN



## 5 LUMEN CVC SETS AND KITS

### “CS” FIVE LUMEN CVC SET

#### SET COMPONENTS

- polyurethane AGB catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-15855	8.5 Fr x 8"	K	.032" J/SST		X	X	5

### “CS” FIVE LUMEN ARROWg<sup>ard</sup> BLUE (A+g) CVC SET

#### SET COMPONENTS

- polyurethane AGB catheter with Blue FlexTip
- guide wire w/ Advancer
- insertion needle 18 Ga x 6.35 cm
- dilator
- w/ or w/o Raulerson syringe 5 ml

ARROW

REF.	CATHETER	LUMEN	SWG	RAULERSON	DILATOR	CE	QTY
CS-25855	8.5 Fr x 8"	K	.032" J/SST		X	X	5
CS-25855-E	8.5 Fr x 8"	K	.032" J/SST	X	X	X	5



# KITS CONFIGURED TO MEET YOUR PROCEDURAL NEEDS

## THEY SAVE TIME, INCREASE EFFICIENCY AND REDUCE WASTE

Individual solutions for optimum patient care, combined with ease of handling. Our European ARROW Select Kits (EASK) are the perfect way to ensure compliance combined with safety, easy handling and cost-effectiveness, and to supply you with your ultimate mix.

EASK simplifies purchasing and reduces storage space requirements for a certain procedure to a single SKU. In addition, there is no more need to stock extra components often not included in a standard kit.

## A PROGRAMME DESIGNED TO SUPPORT YOUR NEEDS

Teleflex has formed a multifunctional team of EASK specialists, who focus on processing and expediting your EASK orders.

The EASK team can also provide a no-obligation analysis to identify the best configuration of products for your customised set. We want to meet your particular requirements.

### References:

- 1 Veenstra, D.L., Saint, S., Saha, S., Lumley, T., Sullivan, S.D. "Cost-Effectiveness of Antiseptic-Impregnated Central Venous Catheters for the Prevention of Catheter-Related Bloodstream Infection," *Journal of the American Medical Association*, January 20, 1999, Vol. 281, Issue 3, pp. 261-267.
- 2 Maki, D.G., Stolz, S.M., Wheeler, S., Mermel, L.A. "Prevention of Central Venous Catheter-Related Bloodstream Infection With an Antiseptic-Impregnated Catheter: A Randomized, Controlled Trial," *Annals of Internal Medicine*, August 15, 1997, Vol. 127, Issue 4, pp. 257-266.
- 3 References include but are not limited to Maki 1997, Sampath 1995 and Bach 1994.
- 4 O'Grady et al. "Guidelines for the Prevention of Intravascular Catheter-Related Infections," *Infection Control Hospital Epidemiology*, Dec. 2002, Vol. 23, Issue 12, pp. 759-69.

Full bibliography available upon request.



Our European ARROW Select Kits (EASK) programme enables you to order special catheter kits configured to meet your requirements.

## EASK EUROPEAN ARROW SELECT KITS

- economical & ecological – save time and costs, avoid waste
- geared to your internal operational structures and preferences
- improve safety for patient and caregiver

Teleflex Incorporated (NYSE: TFX) is a diversified global company, distinguished by a significant presence in healthcare, with niche businesses that also serve the aerospace and commercial markets.

Teleflex Medical is committed to partnering with healthcare providers in anaesthesiology, critical care, urology and surgery to provide solutions that help reduce infections and improve patient and provider safety. The company also produces surgical instruments and devices, cardiac devices and other specialty products for device manufacturers.

The Teleflex Medical family of brands includes ARROW<sup>®</sup>, BEERE<sup>®</sup>, DEKNATEL<sup>®</sup>, GIBECK<sup>®</sup>, HUDSON RCI<sup>®</sup>, KMEDIC<sup>®</sup>, PILLING<sup>®</sup>, PLEUR-EVAC<sup>®</sup>, RÜSCH<sup>®</sup>, SHERIDAN<sup>®</sup>, SMD<sup>™</sup>, TAUT<sup>®</sup>, TFX OEM<sup>®</sup> and WECK<sup>®</sup>.

Teleflex Medical global operations: Austria, Belgium, Canada, China, Czech Republic, France, Germany, Greece, Hungary, India, Ireland, Italy, Japan, Mexico, Netherlands, Portugal, Singapore, Slovak Republic, South Africa, Spain, Switzerland, United Kingdom, Uruguay and USA.

#### **YOUR CONTACTS FOR EUROPE, MIDDLE EAST AND AFRICA (EMEA):**

##### **TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND**

Teleflex Medical Europe Ltd., IDA Business Park, Athlone, Co Westmeath  
Phone +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73  
orders.intl@teleflexmedical.com

**AUSTRIA** +43 (0)1 402 47 72  
**BELGIUM** +32 (0)2 333 24 60  
**CZECH REPUBLIC** +420 (0)495 759 111  
**FRANCE** +33 (0)5 62 18 79 40  
**GERMANY** +49 (0)7151 406 0  
**GREECE** +30 210 67 77 717  
**HUNGARY** +36 (0)1 475 13 60  
**ITALY** +39 0362 58 911  
**NETHERLANDS** +31 (0)88 00 215 00  
**PORTUGAL** +351 22 541 90 85  
**SLOVAK REPUBLIC** +421 (0)3377 254 28  
**SOUTH AFRICA** +27 (0)11 807 4887  
**SPAIN** +34 918 300 451  
**SWITZERLAND** +41 (0)31 818 40 90  
**UNITED KINGDOM** +44 (0)1494 53 27 61

For detailed information see [www.teleflexmedical.com/contact](http://www.teleflexmedical.com/contact)

The products in this catalogue are only available for EMEA (Europe, Middle East, Africa). For further information contact your local representative. All data current at time of printing (09/2010). Subject to technical changes without further notice.

94 10 78 - 00 00 21 · REV A · MC / WM · 08 10 05