


Arrow

# Integrated Arterial Catheter

## The Original All-In-One

With its over-needle design, the Arrow Integrated Arterial Catheter allows for quick and easy insertions with high success rates. Often referred to as "the dart" because of its design, it was the first integrated arterial catheter offered. Today, it continues to deliver on its leading-edge innovation.

- **Saves Time:** Integrated all-in-one design
- **Delivers Safety:** Designed to minimise blood exposure
- **Reduces Complications:** Less risk of dislodgement, transfixing, and losing cannulation



# Ordering Information

## Arrow Integrated Arterial Catheter

ITEM NUMBER	CATHETER			GUIDEWIRE			INTRODUCER NEEDLE		REFERENCE	TYPE	QTY/ CASE
	DIA (G)	LENGTH CM (IN)	MATERIAL	DIA MM (IN)	LENGTH (CM)	TIP CONFIG.	DIA (G)	LENGTH CM (IN)	NOTES		
RA-04122	22	3.49 (1.37")	PUR	0.38 (.015")	12.4	Straight	23*	8.22 (3.24")	Integral Guidewire (2-piece)	Set	50
RA-04120	20	3.81 (1.50")	PUR	0.46 (.018")	13.7	Straight	22*	9.17 (3.61")	Integral Guidewire	Set	50
RA-04020	20	4.45 (1.75")	PUR	0.46 (.018")	13.7	Straight	22*	9.66 (3.80")	Integral Guidewire	Set	50

\*TW: Thin Wall