


Arrow

Seldinger Arterial Catheter

Advanced design for advanced performance.

The Arrow Seldinger Arterial Catheter is designed to improve patient safety by eliminating confusion of catheter identification, reduce complications associated with insertion technique, and provide optimal performance leading to effective treatment for patients. Available in multiple sizes to provide customisation to each patient, each insertion site, and each application.

- **Designed for Patient Safety:** Minimises trauma to vessel wall and eliminates confusion of catheter identification
- **Reduces Complications:** Designed for less vessel movement, irritation, infiltration, and accidental disconnections
- **Saves Time and Costs:** Eliminates finding and purchasing an add-on extension set
- **Advanced Performance:** Innovative catheter material demonstrated high resistance to collapse, low thrombus accumulation and occlusion¹



Ordering Information

Arrow Seldinger Arterial Catheter Sets'

ITEM NUMBER	CATHETER			GUIDEWIRE					INTRODUCER NEEDLE		QTY/CASE
	DIA (G)	LENGTH CM (IN)	PUR	DIA MM (IN)	LENGTH (CM)	TIP CONFIGURATION	EXIT MARKING	REFERENCE	DIA (G)	LENGTH CM (IN)	
SAC-00520	20	5 (1.97")	PEBA	0.53 (.021")	35	Straight/Straight	•		20**	4 (1.57")	10
SAC-00820	20	8 (3.15")	PEBA	0.53 (.021")	35	Straight/Straight	•		20**	4 (1.57")	10
SAC-01220	20	12 (4.72")	PEBA	0.53 (.021")	35	Straight/Straight	•		20**	7 (2.76")	10
SAC-00818	18	8 (3.15")	PEBA	0.64 (.025")	33	Straight/Straight	•		18**	5 (1.97")	10
SAC-01218	18	12 (4.72")	PEBA	0.64 (.025")	33	Straight/Straight	•		18**	5 (1.97")	10

*Latex-free, PVC-free, DEHP-free **XTW: Extra Thin Wall

Rx only

References: 1. Teleflex. Seldinger Arterial Catheter Benchtop Studies. Data on File, D011453: 2014. 2. Teleflex. Testing of Introducer Needles. Data on File, AIMH-001296: 2014.

Teleflex, the Teleflex logo and Arrow are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are trademarks of their respective owners. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current as of 10/2020. © 2020 Teleflex Incorporated. All rights reserved. MCI-100033-EN-IN Rev 0 10 20